





Intubation Aids

Stylets / Introducers / Tube exchanger

S-Guide

Malleable intubating guide for difficult airway management

Features



Soft and Atraumatic Tip

- Preformed
- Coloured:
 - For better visualisation within the glottic structures
 - For easy size recognition

Oxygenation* (11 Fr and 15 Fr)

- Three openings at the distal end of the S-Guide to ensure appropriate oxygen flow
- Supplied with two types of O₂ Connectors

Malleability

- Malleability and rigidity, to preform and manoeuvre the tracheal tube
- Black marking on the sizes 11 Fr and 15 Fr shows ideal position to preform the tracheal tube

Order
information

S-Guide / Two types of O₂ Connectors included (11 Fr and 15 Fr), for single use, sterile

| Size | Length | Centimetre Markings | For Tracheal Tube | REF | Box |
|-------|--------|---------------------|-------------------|--------------------|-----|
| 8 Fr | 40 cm | 5 to 30 cm | ≥ 3.0 mm I.D. | 33-92-400-1 | 5 |
| 11 Fr | 65 cm | 5 to 40 cm | ≥ 4.5 mm I.D. | 33-91-650-1 | 5 |
| 15 Fr | 65 cm | 5 to 40 cm | ≥ 6.0 mm I.D. | 33-90-650-1 | 5 |

*Oxygenation is possible with the sizes 11 Fr and 15 Fr. Oxygenation is not possible with the size 8 Fr.

Additional
information

Application video



Brochure



Recommendation of use



Annotated bibliography

/ i-Bougie

The versatile Introducer to facilitate difficult intubation

Features




Atraumatic Distal Tip

- Rounded and angled design with lateral openings




Oxygenation

- Allows oxygen flow through the lateral openings at the distal tip
- Supplied with two types of O₂ Connectors



Shapeability

- Material allows shaping of the i-Bougie
- Centimetre markings from 5 to 40 cm
- Low friction surface for easy railroading of the tracheal tube



Order
information

i-Bougie / Two types of O₂ Connectors included, for single use, sterile

| Size | Length | Centimetre Markings | For Tracheal Tube | REF | Box |
|-------|--------|---------------------|-------------------|-------------|-----|
| 14 Fr | 70 cm | 5 to 40 cm | ≥ 6.0 mm I.D. | 33-70-700-1 | 5 |

Additional
information



Application video



Brochure



Annotated bibliography

X-Changer

Tube exchanger to maintain the airway during tracheal tube exchange




Features



Soft and Atraumatic Tip

- Blunt tip is flexible and atraumatic
- Colour coded for easy size recognition



Oxygenation

- Hollow catheter with lateral openings to deliver oxygen
- Supplied with two types of O₂ Connectors



Centimetre Markings

- For safe determination of the insertion depth
- Allow proper alignment with the tracheal tube

3 mm



Length 80 cm

VBM

CE 0123

Order
information**X-Changer** / Two types of O₂ Connectors included, for single use, sterile

| Size | Length | Centimetre Markings | For Tracheal Tube | REF | Box |
|------|--------|---------------------|-------------------|--------------------|-----|
| 11Fr | 80 cm | 5 to 40 cm | ≥ 4.0 mm I.D. | 33-42-800-1 | 5 |
| 14Fr | 80 cm | 5 to 40 cm | ≥ 5.0 mm I.D. | 33-41-800-1 | 5 |
| 19Fr | 80 cm | 5 to 40 cm | ≥ 7.0 mm I.D. | 33-40-800-1 | 5 |

Additional
information

Application video





Brochure



Annotated bibliography

/ Pocket Introducer

The smallest packaged introducer to facilitate difficult intubation



Features

Tip

- Preformed and rounded olive tip is atraumatic
- Distal end is shapeable

Click Connections

- Folded to 20 cm, unfolds to 65 cm with only 3 clicks
- The click connections are permanent and cannot be detached again
- Centimetre markings from 10 to 40 cm

Smallest Pack Size

- Ideal space solution for emergency bags, helicopters, lab coat pockets, crash carts, etc.
- Ready for use within seconds



Order information

Pocket Introducer / For single use, sterile

| Size | Length | Centimetre Markings | For Tracheal Tube | REF | Box |
|-------|--------|---------------------|-------------------|--------------------|-----|
| 15 Fr | 65 cm | 10 to 40 cm | ≥ 6.0 mm I.D. | 33-80-650-1 | 5 |

Pocket Introducer / In robust plastic cartridge, for single use, sterile

| Size | Length | Centimetre Markings | For Tracheal Tube | REF | Box |
|-------|--------|---------------------|-------------------|--------------------|-----|
| 15 Fr | 65 cm | 10 to 40 cm | ≥ 6.0 mm I.D. | 33-81-650-1 | 1 |

Additional information



Application video




Annotated bibliography



Plastic cartridge

/ Stylet

Intubation stylet enables preforming of the tracheal tube



Features

Tip

- Soft, atraumatic tip

Versions

- Without silicone connector
- With moveable silicone connector, which allows fixation of the tracheal tube

Order
information**Stylet / Reusable**

| Size | Length | For Tracheal Tube | Without Silicone Connector | Box | With Moveable Silicone Connector (15 mm I.D.) | Box |
|--------|--------|-------------------|----------------------------|-----|---|-----|
| Small | 34 cm | 2.5–3.5 mm I.D. | REF 90-10-111 | 25 | REF 90-10-001 | 10 |
| Medium | 38 cm | 4.0–4.5 mm I.D. | REF 90-10-222 | 25 | REF 90-10-002 | 10 |
| Large | 45 cm | 5.0–6.5 mm I.D. | REF 90-10-333 | 25 | REF 90-10-003 | 10 |
| Large | 45 cm | 7.0–11.0 mm I.D. | REF 90-10-444 | 25 | REF 90-10-004 | 10 |

Additional
information

Annotated bibliography

Additional
information



Find our full range of products on www.vbm-medical.de

Here you will find our media centre with additional information such as application videos, brochures, flyer etc.

The medical devices in this advertising material are manufactured without the use of natural rubber latex, unless otherwise specified.
The medical devices in this advertising material do not contain phthalates which require labelling according to CLP Regulation (EC) 1272/2008.

VBM *Medizintechnik GmbH*

Einsteinstrasse 1 / 72172 Sulz a. N. / Germany

Tel.: +49 7454 9596-0 / Fax: +49 7454 9596-33 / e-mail: sales@vbm-medical.de / www.vbm-medical.de

Connect with us:



CE 0123
FLY013_5.0_EN