


## Accessories for Anaesthesia & Intensive Care



Tube Fixations / Respiration Tubing / Respiration Accessories  
Catheter Mount Tubing / Connectors / Smooth Silicone Tubing  
Pressure Infusion / Support Arms / Tube Connectors



Find our full range of products on [www.vbm-medical.de](http://www.vbm-medical.de)

Here you will find our media centre with additional information such as application videos, brochures, flyers etc.

**Disclaimer:**

This advertising material includes information of the complete VBM medical device product range.

Whether a medical device is available / approved in your region must be inquired. Please contact your VBM customer service.

The medical devices in this catalogue are manufactured without the use of natural rubber latex, unless otherwise specified.

The medical devices in this catalogue do not contain phthalates which require labelling according to CLP Regulation (EC) 1272/2008.

---

	<a href="#">Tube Fixations</a>
Tracheostomy Necktapes	<b>Page 4 – 5</b>
Biteblock	
Endotape	
Tube Holder	
	<a href="#">Respiration Tubing</a>
Respiration Tubing for Adults	<b>Page 6</b>
Respiration Tubing for Children	
Anaesthesia System for Children	
	<a href="#">Respiration Accessories</a>
Rebreathing Bags	<b>Page 7</b>
Water Traps	
	<a href="#">Catheter Mount Tubing</a>
Catheter Mount Tubing Reusable	<b>Page 8 – 9</b>
Connecting Tubing	
Catheter Mount Tubing for Single Use	
	<a href="#">Connectors</a>
Straight Connectors	<b>Page 10 – 14</b>
Right Angled Connectors	
Tube Inhaler	
Elbows	
Y-Connectors	
Paediatric Y-Connectors	
T-Connectors	
Paediatric Connecting Tubing	
Accessories for Connectors	
	<a href="#">Smooth Silicone Tube</a>
Smooth Silicone Tube	<b>Page 15</b>
	<a href="#">Pressure Infusion</a>
Classic Infusor	<b>Page 16 – 19</b>
Clear Infusor	
Dispo Infusor	
Accessories for Pressure Infusion Cuffs	
	<a href="#">Support Arms</a>
Support Arm with Flexible Rod	<b>Page 20 – 21</b>
Support Arm with 3-Joint-System	
Tube-Tray	
Pressure Transducer Holders	
	<a href="#">Tube Connectors</a>
Positive Locking Connectors	<b>Page 22</b>
Luer Lock Connectors	

# Tracheostomy Necktapes

Tube fixations for tracheostomised patients

Order information

**Tracheostomy Necktape** / For single use, colour: blue

Patient	Length	REF	Box
Infant / Child	355 mm	40-0005-011	10
Child / Adult	475 mm	40-0005-033	10
Adult, extra large	525 mm	40-0005-044	10



**Tracheostomy Necktape** / For single use, colour: beige

Patient	Length	REF	Box
Infant / Child	340 mm	40-0008-011	10
Child / Adult	490 mm	40-0008-033	10
Adult, extra large	590 mm	40-0008-044	10



**Tracheostomy Necktape, Separable** / For single use, colour: beige

Patient	Length	REF	Box
Child / Adult	400 mm	40-0066-033	10



# Biteblock


Oral tube fixation for tracheal tube, Laryngeal Tube and laryngeal mask

Tube clips into biteblock and is therefore protected from occlusion.

Order information

**Biteblock** / For single use, sterile

Size	Tracheal Tube	Laryngeal Tube	Laryngeal Mask	Biteblock	Biteblock with Elastic Necktape	Box
4	6.5-8.0 mm I.D.	-	2, 2.5	REF 40-60-040-1	REF 40-66-040-1	20
5	8.5-11.0 mm I.D.	2, 2.5	-	REF 40-60-050-1	REF 40-66-050-1	20
6	-	3, 4, 5	3, 4, 5	REF 40-60-060-1	REF 40-66-060-1	20



# Endotape

Oral/Nasal tracheal tube fixation

Order information

**Endotape** / With elastic, silicone coated ties, for single use, colour: beige

Patient	Length	REF	Box
Child / Adult	540 mm	40-8018-044*	10

\* Contains or presence of natural rubber latex



Additional information



Application video

# Tube Holder

Oral/Nasal tracheal tube fixation

Order information

**Tube Holder** / For single use, colour: beige

Patient	Length	REF	Box
Child / Adult	540 mm	40-8010-033	10
Adult	650 mm	40-8010-044	10



Additional information




Application video

# Respiration Tubing 19 mm I.D. for Adults

Reusable

Order information

**Vulcanised 22 mm I.D. Cuffs at Both Ends /**



Length	Colour: green	Box
30 cm	REF 66-22-030	1
40 cm	REF 66-22-040	1
60 cm	REF 66-22-060	1
80 cm	REF 66-22-080	1

Length	Colour: green	Box
90 cm	REF 66-22-090	1
110 cm	REF 66-22-110	1
130 cm	REF 66-22-130	1
150 cm	REF 66-22-150	1


Length	Colour: green	Box
180 cm	REF 66-22-180	1
200 cm	REF 66-22-200	1

# Respiration Tubing 9 mm I.D. for Children

Reusable

Order information

**Vulcanised 10 mm I.D. Cuffs at Both Ends /**




Length	Colour: green	Box
25 cm	REF 66-10-025	1
45 cm	REF 66-10-045	1
60 cm	REF 66-10-060	1

Length	Colour: green	Box
85 cm	REF 66-10-085	1
110 cm	REF 66-10-110	1
130 cm	REF 66-10-130	1

Length	Colour: green	Box
150 cm	REF 66-10-150	1
200 cm	REF 66-10-200	1

**Vulcanised 10 mm / 22 mm I.D. Cuffs /**




Length	Colour: green	Box
25 cm	REF 66-11-025	1
45 cm	REF 66-11-045	1
60 cm	REF 66-11-060	1

Length	Colour: green	Box
85 cm	REF 66-11-085	1
110 cm	REF 66-11-110	1
130 cm	REF 66-11-130	1

Length	Colour: green	Box
150 cm	REF 66-11-150	1
200 cm	REF 66-11-200	1

**Vulcanised 22 mm I.D. Cuffs at Both Ends /**



Length	Colour: green	Box
25 cm	REF 66-12-025	1
60 cm	REF 66-12-060	1
85 cm	REF 66-12-085	1

Length	Colour: green	Box
110 cm	REF 66-12-110	1
130 cm	REF 66-12-130	1

Length	Colour: green	Box
150 cm	REF 66-12-150	1
200 cm	REF 66-12-200	1

# Anaesthesia System for Children

Reusable

Order information

**Anaesthesia System /**

- 2 x 9 mm I.D. Silicone Respiration Tubing, 22 mm / 10 mm I.D. Cuff, 110 cm long
- 1 x 9 mm I.D. Silicone Respiration Tubing, 22 mm / 22 mm I.D. Cuff, 110 cm long
- 1 x Right Angled Paediatric Y-Connector, standard, swivelling, 15 mm I.D. / 22 mm O.D. - 10 mm O.D. - 10 mm O.D.
- 1 x Straight Connector
- 1 x Silicone Rebreathing Bag 500 ml

REF	Box
62-11-200	1



# Rebreathing Bags

## Reusable

Order  
information

### Standard, 22 mm I.D. Cuff and Loop /

Size	REF	Box
500 ml	66-00-011	1
1000 ml	66-00-015	1
1500 ml	66-00-012	1
2300 ml	66-00-013	1
3000 ml	66-00-014	1

22 mm I.D.



### Open Tube Connection and Cap, 22 mm I.D. Cuff and Loop /

Size	REF	Box
500 ml	66-00-041	1
1000 ml	66-00-045	1
1500 ml	66-00-042	1
2300 ml	66-00-043	1
3000 ml	66-00-044	1

22 mm I.D.



Tube connection and cap

### 5 mm Hole, 22 mm I.D. Cuff and Loop /

Size	REF	Box
500 ml	66-00-441	1
1000 ml	66-00-445	1
1500 ml	66-00-442	1
2300 ml	66-00-443	1
3000 ml	66-00-444	1

22 mm I.D.



5 mm hole

# Water Traps

## Reusable

Order  
information

**Water Trap / Adult, 15 mm I.D. / 22 mm O.D., container 55 ml**

REF	Box
66-40-500	6

15 mm I.D.  
22 mm O.D.



55 ml

15 mm I.D.  
22 mm O.D.

**Water Trap / Child, with tube connection 10 mm O.D., container 55 ml**

REF	Box
66-40-510	6

10 mm O.D.



55 ml

10 mm O.D.

Accessories  
for Respiration  
Tubing

**Hook / For respiration tubing 9 mm I.D. / 19 mm I.D., reusable**

REF	Box
66-00-002	5



# Catheter Mount Tubing

Reusable

Order information

## Straight Connector /

Length: 13 cm	Length: 20 cm	Box
REF 65-50-130	REF 65-50-200	10


Length: 13 cm	Length: 20 cm	Box
REF 65-55-130	REF 65-55-200	10



## Right Angled Swivel Connector /

Length: 13 cm	Length: 20 cm	Box
REF 65-60-130	REF 65-60-200	10

Length: 13 cm	Length: 20 cm	Box
REF 65-80-130	REF 65-80-200	10



## Right Angled Swivel Connector, Double Swiveling, with Bronchoscope/Suction Cap /

Length: 13 cm	Length: 20 cm	Box
REF 65-68-130	REF 65-68-200	10

Length: 13 cm	Length: 20 cm	Box
REF 65-88-130	REF 65-88-200	10



# Connecting Tubing

Reusable

Order information

Length: 13 cm	Length: 20 cm	Box
REF 65-00-130	REF 65-00-200	10





# Catheter Mount Tubing

For single use, sterile

Order  
information

## Straight Connector /

Length: 15 cm Box

REF 60-50-150-1 25

Length: 15 cm Box

REF 60-55-150-1 20

## Straight Swivel Connector /

Length: 15 cm Box

REF 60-52-150-1 25

## Right Angled Swivel Connector /

Length: 15 cm Box

REF 60-60-150-1 25

Length: 15 cm Box

REF 60-80-150-1 20

## Right Angled Swivel Connector with Broncho- scope/Suction Cap /

Length: 15 cm Box

REF 60-65-150-1 25

Length: 15 cm Box

REF 60-85-150-1 20

## Right Angled Swivel Connector, Double Swivel- ing, with Bronchoscope/Suction Cap /


Length: 15 cm Box

REF 60-68-150-1 25

Length: 15 cm Box


REF 60-88-150-1 20

15 mm I.D.  
22 mm O.D.




15 mm O.D.

15 mm I.D.  
22 mm O.D.




22 mm I.D.

15 mm I.D.  
22 mm O.D.




15 mm O.D.

15 mm I.D.  
22 mm O.D.




15 mm O.D.

15 mm I.D.  
22 mm O.D.




22 mm I.D.

15 mm I.D.  
22 mm O.D.




15 mm O.D.

15 mm I.D.  
22 mm O.D.




22 mm I.D.

15 mm I.D.  
22 mm O.D.



15 mm O.D.

15 mm I.D.  
22 mm O.D.



22 mm I.D.

# Straight Connectors

Reusable

Order information

REF	Box
60-30-002	10



REF	Box
60-30-003	10



REF	Box
60-30-005	10



REF	Box
60-40-002	10



REF	Box
60-40-003	10



REF	Box
60-40-005	10



REF	Box
60-12-110	10




REF	Box
60-12-400	10



REF	Box
60-12-200	10



REF	Box
60-12-201	10



REF	Box
60-50-002	10



REF	Box
60-50-005	10




REF	Box
60-50-003	10




REF	Box
60-12-300	10



REF	Box
60-12-601	10



REF	Box
60-71-002	10



# Right Angled Connectors

Reusable

Order  
information

REF	Box
60-65-002	10



REF	Box
60-65-004	10



# Right Angled Swivel Connectors

Reusable

Order  
information


REF	Box
60-60-002	10




REF	Box
60-60-004	10




REF	Box
60-60-010	10




REF	Box
60-60-005	10



REF	Box
60-60-205	10



REF	Box
60-60-305	10



# Right Angled Swivel Connectors


For single use, sterile

Order  
information


REF	Box
60-60-222-1	10



REF	Box
60-60-555-1	10



REF	Box
60-60-111-1	10



# Tube Inhaler

Reusable


For application of MDI canisters on intubated patients.

Order information


REF	Box
60-60-009	10



REF	Box
60-65-009	10



REF	Box
60-60-209	10



REF	Box
60-67-009	10



# Elbows

Reusable

Order information


Reusable	Single use	Box
REF 60-12-000	REF 60-13-000	10



REF	Box
60-12-001	10



REF	Box
60-12-006	10



## Straight Y-Connectors

Reusable

Order information


REF	Box
60-10-100	10




REF	Box
60-11-200	10



REF	Box
60-11-203	10



REF	Box
60-11-204	10



REF	Box
60-11-400	10



REF	Box
60-11-420	10



## Right Angled Y-Connector

Reusable

Order information

REF	Box
60-11-300	10



## Straight Paediatric Y-Connectors


Reusable

Order information

REF	Box
60-20-500	10



REF	Box
60-20-501	10



## Right Angled Paediatric Y-Connectors


Reusable

Order information

REF	Box
60-20-100	10



REF	Box
60-20-101	10



# T-Connectors

Reusable

Order information

REF	Box
60-11-500	10




REF	Box
60-11-503	10



REF	Box
60-11-506	10



REF	Box
60-11-510	10



# Paediatric Connecting Tubing

Reusable

Order information

Length: 9.0 cm	Box
REF 65-00-090	10



Length: 10.5 cm	Box
REF 65-00-105	10




# Accessories for Connectors

Order information

**Meconium Aspirator** / For single use

REF	Box
49-30-000	10



- 15 mm I.D. connection allows direct connection to the tracheal tube for tracheal suctioning of aspirated meconium of the newborn
- Funnel shaped connection for a standard suction catheter

Order information

**Uncoupling Wedge** / Reusable

REF	Box
60-60-999	10



**Sampling Line** / For single use, non-sterile

Version	Length	REF	Box
1.5 mm I.D. / 3.0 mm O.D. with two male Luer Lock	150 cm	61-99-150	10
1.5 mm I.D. / 3.0 mm O.D. with two male Luer Lock	200 cm	61-99-200	10
1.5 mm I.D. / 3.0 mm O.D. with two male Luer Lock	300 cm	61-99-300	10




# Smooth Silicone Tube

Order  
information

Smooth Silicone Tube / Colour: transparent

I.D. (mm)	C (mm)	O.D. (mm)	REF	Box
1.0	0.5	2.0	68-010-05	25 m
1.0	1.0	3.0	68-010-10	25 m
1.5	0.75	3.0	68-015-075	25 m
2.0	1.0	4.0	68-020-10	25 m
3.0	1.0	5.0	68-030-10	25 m
3.0	1.5	6.0	68-030-15	25 m
4.0	1.0	6.0	68-040-10	25 m
4.0	1.5	7.0	68-040-15	25 m
4.0	2.0	8.0	68-040-20	25 m
5.0	1.0	7.0	68-050-10	25 m
5.0	1.5	8.0	68-050-15	25 m
5.0	2.0	9.0	68-050-20	25 m
6.0	1.0	8.0	68-060-10	25 m
6.0	1.5	9.0	68-060-15	25 m
6.0	2.0	10.0	68-060-20	25 m
6.0	3.0	12.0	68-060-30	25 m
7.0	1.5	10.0	68-070-15	25 m
7.0	2.0	11.0	68-070-20	25 m
7.0	3.0	13.0	68-070-30	25 m
8.0	2.0	12.0	68-080-20	25 m
8.0	3.0	14.0	68-080-30	25 m
9.0	2.0	13.0	68-090-20	25 m
10.0	2.0	14.0	68-100-20	25 m
10.0	3.0	16.0	68-100-30	25 m
12.0	3.0	18.0	68-120-30	25 m
12.0	4.0	20.0	68-120-40	25 m



# Classic Infusor

For pressure infusion, irrigation, transfusion and invasive blood pressure monitoring



## 500 ml and 1000 ml Version

With double hook, for fixation of different fluid bag configurations.

- Rapid administration of an infusion solution for compensation of fluid/ volume deficiency
- Pressure monitoring and irrigation of the tubing system during invasive, arterial blood pressure monitoring
- Transfusion

## 3000 / 5000 ml Version

With zipper, allows insertion of different fluid bag configurations.

- Administration of large fluid volumes for irrigation of the operation field

The polyamide material of the Classic Infusor enables easy cleaning & disinfection. Transparent netting for a clear view on the fluid bag and the fluid level. With UDI code for product traceability.



### Features




#### Hand Inflator with Three-Way Stopcock

- Ergonomic hand inflator for inflation
- Three-way stopcock for quick deflation during change of the fluid bag



#### Manometer

- Pressure range from 0 - 300 mmHg in increments of 10 mmHg
- Robust and durable for accurate pressure monitoring



#### Colour Coded Pressure Indicator

- 360° pressure indicator facilitates pressure monitoring
- Pressure relief valve to prevent over inflation

### Order information

#### Classic Infusor / Reusable

Version	500 ml	1000 ml	3000 / 5000 ml	Box
With hand inflator, three-way stopcock and colour coded pressure indicator	REF 58-02-050	REF 58-02-100	-	1
With hand inflator, three-way stopcock and manometer	REF 58-08-050	REF 58-08-100	REF 58-08-300	1
With Luer Lock connection	-	-	REF 58-01-300	1

### Additional information



Application video



Brochure



# / Clear Infusor

For pressure infusion, irrigation, transfusion and invasive blood pressure monitoring



## 500 ml and 1000 ml Version

- Rapid administration of an infusion solution for compensation of fluid/ volume deficiency
- Pressure monitoring and irrigation of the tubing system during invasive, arterial blood pressure monitoring
- Transfusion

## 3000 ml Version

- Administration of large fluid volumes for irrigation of the operation field

Wrap-around design allows fixation of most fluid bag configurations and even pressure distribution. The smooth surface of the PU material enables easy cleaning & disinfection. The transparent material offers a clear 360° view on the fluid bag and fluid level. With UDI code for product traceability.



### Features




#### Hand Inflator With Three-Way Stopcock

- Ergonomic hand inflator for inflation
- Three-way stopcock for quick deflation during change of the fluid bag



#### Manometer

- Pressure range from 0-300 mmHg in increments of 10 mmHg
- Robust and durable for accurate pressure monitoring



#### Colour Coded Pressure Indicator

- 360° pressure indicator facilitates pressure monitoring
- Pressure relief valve to prevent over inflation

### Order information

#### Clear Infusor / Reusable

Version	500 ml	1000 ml	3000 ml	Box
With hand inflator, three-way stopcock and colour coded pressure indicator	REF 57-02-050	REF 57-02-100	-	1
With hand inflator, three-way stopcock and manometer	REF 57-08-050	REF 57-08-100	REF 57-08-300	1
With Luer Lock connection	-	-	REF 57-01-300	1

### Additional information



Application video



Brochure

# / Dispo Infusor

For pressure infusion, irrigation, transfusion and invasive blood pressure monitoring



## 500 ml and 1000 ml Version

With hook fastener, for adjustable fixation of different fluid bag configurations.

- Rapid administration of an infusion solution for compensation of fluid/ volume deficiency
- Pressure monitoring and irrigation of the tubing system during invasive, arterial blood pressure monitoring
- Transfusion

## 3000 ml Version

With hook, for easy fixation of most fluid bag configurations.

- Administration of large fluid volumes for irrigation of the operation field

Transparent netting for a clear view on the fluid bag and the fluid level.




## Features



### Hand Inflator with Three-Way Stopcock

- Ergonomic hand inflator for inflation
- Three-way stopcock for quick deflation during change of the fluid bag



### Colour Coded Pressure Indicator

- 360° pressure indicator facilitates pressure monitoring
- Pressure relief valve to prevent over inflation

## Order information

**Dispo Infusor** / For single patient use

Version	500 ml	1000 ml	3000 ml	Box
With hand inflator, three-way stopcock and colour coded pressure indicator	REF 56-02-050	REF 56-02-100	REF 56-02-300	5
With Luer Lock connection	-	-	REF 56-01-300	5

## Additional information



Application video



Brochure

# Accessories for Pressure Infusion Cuffs

Order information

## Hand Inflator / Reusable

With Three-Way Stopcock and Colour Coded Pressure Indicator	With Three-Way Stopcock	Box
---	-------------------------	-----

REF 51-02-666	REF 51-03-666	1
---------------	---------------	---



With three-way stopcock and colour coded pressure indicator



With three-way stopcock

## Manometer 300 mmHg / Pressure proved, for fixation to pressure infusion cuff, reusable

REF	Box
-----	-----

20-17-301	1
-----------	---



## Support Arm with Flexible Rod

Adjustable in height, length of Support Arm: 38 cm

Order  
information

	Tube Holder Plate	Swivelling Tube Holder Plate	Tube Holder Plate with Safety Lock	Twin Tube Holder	Box
Rail Clamp	REF 73-512-04	REF 73-512-14	REF 73-512-44	REF 73-512-23	1
Universal Clamp	REF 73-532-04	REF 73-532-14	REF 73-532-44	REF 73-532-23	1



## Support Arm with 3-Joint-System

3 joints can be fixed with only one fixation knob

Order  
information

With Universal Clamp and Tube Holder Plate / Available in two lengths

Length	Tube Holder Plate	Tube Holder Plate with Safety Lock	Box
30 cm	REF 70-131-04	REF 70-131-44	1
50 cm	REF 70-231-04	REF 70-231-44	1



## Support Arm with 3-Joint-System

3 joints can be fixed with only one fixation knob

Order  
information

With Universal Clamp and Clamping Piece for Transducer Plates / Available in two lengths

Length	REF	Box
30 cm	70-131-07	1
50 cm	70-231-07	1



# / Tube-Tray

Tube holders for adult and paediatric respiration tubing



- Various openings for a perfect fit of respiration tubing of different size and material
- Also accepts gas sampling tubes
- Does not contain ferromagnetic materials (except Tube-Tray with rail clamp)
- Versatile placement either behind or lateral to the patients' head
- Relieves stress on patients' face

Order  
information

## Tube-Tray Straight /

REF	Box
77-10-000	1



## Tube-Tray Angled /

REF	Box
77-20-000	1



## Tube-Tray with Rail Clamp /

REF	Box
77-00-020	1



# / Pressure Transducer Holders

For pressure transducer plates and stopcocks

Order  
information

## Stopcock Holder /

REF	Box
78-10-000	1



## Standard Holder /

REF	Box
78-20-000	1



## Universal Holder with One Flexible Joint /

REF	Box
78-30-000	1





# Positive Locking Connectors

Order information

REF	Box
SLZM30	5

REF	Box
SLZM40	5



REF	Box
SLZF30	5

REF	Box
SLZF40	5

# Luer Lock Connectors

Order information

REF	Box
CHLLM05	5


REF	Box
CHLLM06	5

REF	Box
CHLLM30	5

REF	Box
CHLLM40	5

REF	Box
CHLLMB20	5

REF	Box
CHLLMB30	5



REF	Box
CHLLF05	5

REF	Box
CHLLF06	5

REF	Box
CHLLF30	5

REF	Box
CHLLF40	5

REF	Box
CHLLFB20	5

REF	Box
CHLLFB30	5

<b>A</b>		<b>P</b>	
Accessories for Connectors	14	Paediatric Connecting Tubing	14
Accessories for Pressure Infusion Cuffs	19	Paediatric Y-Connectors	13
Accessories for Respiration Tubing	7	Pressure Infusion	16–19
Anaesthesia System for Children	6	Pressure Infusion Cuffs	16–18
Aspirator	14	Positive Locking Connectors	22
		Pressure Transducer Holders	21
<b>B</b>		<b>R</b>	
Biteblock	4	Rebreathing Bags	7
		Respiration Accessories	7
<b>C</b>		Respiration Tubing for Adults	6
Catheter Mount Tubing	8–9	Respiration Tubing for Children	6
Classic Infusor	16	Right Angled Connectors	11
Clear Infusor	17		
Connectors	10–14	<b>S</b>	
Connectors Accessories	14	Sampling Line	14
Connectors Luer Lock	22	Smooth Silicone Tubing	15
Connectors Straight	10	Straight Connectors	10
Connectors Right Angled	11	Support Arms	20–21
Connectors-T	14	Support Arm with 3-Joint-System	20
Connecting Tubing	8	Support Arm with Flexible Rod	20
Connectors-Y	13		
Connectors-Y for Paediatric	13	<b>T</b>	
		T-Connectors	14
<b>D</b>		Tracheostomy Necktapes	4
Dispo Infusor	18	Tube Connectors	22
		Tube Fixations	4–5
<b>E</b>		Tube Holder	5
Elbows	12	Tube Inhaler	12
Endotape	5	Tube-Tray	21
		<b>U</b>	
<b>H</b>		Uncoupling Wedge	14
Hand Inflator	19		
Hook	7	<b>W</b>	
		Water Traps	7
<b>L</b>			
Luer Lock Connectors	22	<b>Y</b>	
		Y-Connectors	13
<b>M</b>		Y-Connectors for Paediatric	14
Manometer	19		
Meconium Aspirator	14		
		<b>N</b>	
<b>N</b>		Necktapes	4



## **VBM** Medizintechnik GmbH

Einsteinstrasse 1  
72172 Sulz a. N.  
Germany  
Tel.: +49 7454 9596-0  
Fax: +49 7454 9596-33  
e-mail: [sales@vbm-medical.de](mailto:sales@vbm-medical.de)  
[www.vbm-medical.de](http://www.vbm-medical.de)

More product catalogues, brochures, flyers and application videos can be found in the media library on our website.

Connect with us:

