

ORAL AIRWAY

LATEX FREE | SINGLE USE | DISPOSABLE

Latex free oral Guedel airways are designed with built-in bite blocks. Both styles are color coded for accurate and quick size identification. These airways are non-sterile and individually packaged; sets are ideal for crash carts.



COLOR CODED GUEDEL AIRWAY

- Firm, polyethylene plastic
- Built-in bite block
- Individually wrapped

ITEM	DESCRIPTION	COLOR	SIZE	PK
1-1504-40	Newborn	Pink	4 cm / 40 mm	10
1-1504-50	Infant	Blue	5 cm / 50 mm	10
1-1504-60	Small Child	Black	6 cm / 60 mm	10
1-1504-70	Child	White	7 cm / 70 mm	10
1-1504-80	Small Adult	Green	8 cm / 80 mm	10
1-1504-90	Medium Adult	Yellow	9 cm / 90 mm	10
1-1504-99	Large Adult	Red	10 cm / 100 mm	10
1-1504-11	X-Large Adult	Orange	11 cm / 110 mm	10
1-1505-40	Newborn	Pink	4 cm / 40 mm	50
1-1505-50	Infant	Blue	5 cm / 50 mm	50
1-1505-60	Small Child	Black	6 cm / 60 mm	50
1-1505-70	Child	White	7 cm / 70 mm	50
1-1505-80	Small Adult	Green	8 cm / 80 mm	50
1-1505-90	Medium Adult	Yellow	9 cm / 90 mm	50
1-1505-99	Large Adult	Red	10 cm / 100 mm	50
1-1505-11	X-Large Adult	Orange	11 cm / 110 mm	50

SETS

ITEM	SIZES INCLUDED	PK
1-1504-00	40, 50, 60, 70, 80, 90, 100, 110 mm (1 of each size per set/pack)	1

COLOR CODED SOFT PVC GUEDEL AIRWAY

- Soft, clear PVC plastic
- Built-in bite block
- Individually wrapped

ITEM	DESCRIPTION	COLOR	SIZE	PK
1-5010-30	Premature	Clear	3 cm / 30 mm	10
1-5010-20	Newborn	Pink	4 cm / 40 mm	10
1-5010-10	Infant	Blue	5 cm / 50 mm	10
1-5010-01	Small Child	Black	6 cm / 60 mm	10
1-5010-02	Child	White	7 cm / 70 mm	10
1-5010-03	Small Adult	Green	8 cm / 80 mm	10
1-5010-04	Medium Adult	Yellow	9 cm / 90 mm	10
1-5010-05	Large Adult	Red	10 cm / 100 mm	10
1-5010-06	X-Large Adult	Orange	11 cm / 110 mm	10
1-1521-30	Premature	Clear	3 cm / 30 mm	50
1-1521-20	Newborn	Pink	4 cm / 40 mm	50
1-1521-10	Infant	Blue	5 cm / 50 mm	50
1-1521-01	Small Child	Black	6 cm / 60 mm	50
1-1521-02	Child	White	7 cm / 70 mm	50
1-1521-03	Small Adult	Green	8 cm / 80 mm	50
1-1521-04	Medium Adult	Yellow	9 cm / 90 mm	50
1-1521-05	Large Adult	Red	10 cm / 100 mm	50
1-1521-06	X-Large Adult	Orange	11 cm / 110 mm	50

SETS

ITEM	SIZES INCLUDED	PK
1-5010-98	30, 40, 50, 60, 70, 80, 90, 100, 110 mm (1 of each size per set/pack)	1

ORAL AIRWAY

LATEX FREE | SINGLE USE | DISPOSABLE

AirLife™'s latex free oral airways are made of firm plastic with built-in bite blocks. The Berman airways are color coded for accurate and quick size identification. These single use airways are non-sterile; sets are ideal for crash carts.



COLOR CODED BERMAN AIRWAY

- Firm, polyethylene plastic
- Vented septum
- Individually wrapped

ITEM	DESCRIPTION	COLOR	SIZE	PK
1-1508-40	Newborn	Pink	4 cm / 40 mm	10
1-1508-50	Infant	Blue	5 cm / 50 mm	10
1-1508-60	Small Child	Black	6 cm / 60 mm	10
1-1508-70	Child	White	7 cm / 70 mm	10
1-1508-80	Small Adult	Green	8 cm / 80 mm	10
1-1508-90	Medium Adult	Yellow	9 cm / 90 mm	10
1-1508-99	Large Adult	Red	10 cm / 100 mm	10
1-1508-11	X-Large Adult	Orange	11 cm / 110 mm	10

ITEM	DESCRIPTION	COLOR	SIZE	PK
1-1506-40	Newborn	Pink	4 cm / 40 mm	50
1-1506-50	Infant	Blue	5 cm / 50 mm	50
1-1506-60	Small Child	Black	6 cm / 60 mm	50
1-1506-70	Child	White	7 cm / 70 mm	50
1-1506-80	Small Adult	Green	8 cm / 80 mm	50
1-1506-90	Medium Adult	Yellow	9 cm / 90 mm	50
1-1506-99	Large Adult	Red	10 cm / 100 mm	50
1-1506-11	X-Large Adult	Orange	11 cm / 110 mm	50

SETS

ITEM	SIZES INCLUDED	PK
1-1508-00	40, 50, 60, 70, 80, 90, 100, 110 mm (1 of each size per set; 10 sets per pack)	10

WILLIAMS AIRWAY

- Guides fiber optic bronchoscopy
- One piece construction
- Bulk packaged

ITEM	DESCRIPTION	COLOR	SIZE	PK
1-1503-09	Adult Female	Yellow	9 cm	10
1-1503-10	Adult Male	Pink	10 cm	10

