

NOVO CP Cuff Pressure Monitor 30



Cuff pressure manometer

Product properties

- Monitor for measuring and checking the cuff pressure of endotracheal tubes and tracheal cannulae
- Luer-Slip male
- Suitable for MRI
- Single-patient product

Patient target group

- Adolescents
- Adults



REF	Article	Description	SU
N51710030	NOVO CP Cuff Pressure Monitor 30	Cuff pressure manometer	10 pieces



Novo Klinik-Service GmbH
Technologiepark West – Haus 15
Zum Frenser Feld 1
50127 Bergheim • Germany

Rev.01 Novo-update-CPcuff_EN_202405

Tel. +49 2271 99404-0
Fax +49 2271 99404-700
info@novo-med.de
novo-med.de

