


Available in 4 sizes

Autoclavable and single use sterile versions

Powder-free



Esmarch Bandages

For exsanguination of upper and lower extremities

/ Esmarch Bandages

Exsanguinating the extremity before tourniquet cuff inflation facilitates a relatively bloodless field during surgical interventions. The bandages are available in 4 different sizes to ensure easy wrapping according to different patient demands.

Features



Material

- Textured surface for a good grip
- Elastic properties for easy wrapping
- Robust and tear resistant
- Made of synthetic rubber and powder-free (talcum-free) to avoid possible allergic reactions

Order information

Esmarch Bandage / For single use, sterile packed, colour: light blue

Width	Length	REF	Box
6 cm	365 cm	20-21-060-1	10
8 cm	365 cm	20-21-080-1	10
10 cm	365 cm	20-21-100-1	10
15 cm	365 cm	20-21-150-1	10

Esmarch Bandage / Reusable, autoclavable, colour: dark blue

Width	Length	REF	Box
6 cm	500 cm	20-20-060	5
8 cm	500 cm	20-20-080	5
10 cm	500 cm	20-20-100	5
15 cm	500 cm	20-20-150	5

This product is manufactured without the use of natural rubber latex, unless otherwise specified.
This product does not contain phthalates which require labelling according to CLP Regulation (EC) 1272/2008.