

Reduces risk of contamination

Latex free



Dispo Infusor

Pressure infusion cuff, for single patient use

/ Dispo Infusor

500 ml and 1000 ml Version

With hook fastener, for adjustable fixation of different fluid bag configurations.

- Rapid administration of an infusion solution for compensation of fluid/volume deficiency
- Pressure monitoring and irrigation of the tubing system during invasive, arterial blood pressure monitoring
- Transfusion

3000 ml Version

With hook, for easy fixation of most fluid bag configurations.

- Administration of large fluid volumes for irrigation of the operation field


Transparent netting for a clear view on the fluid bag and the fluid level.

Features



Hand Inflator with Three-Way Stopcock

- Ergonomic hand inflator for inflation
- Three-way stopcock for quick deflation during change of the fluid bag



Colour Coded Pressure Indicator

- 360° pressure indicator facilitates pressure monitoring
- Pressure relief valve to prevent over inflation



Order information

Dispo Infusor / For single patient use

Version	500 ml	1000 ml	3000 ml	Box
With hand inflator, three-way stopcock and colour coded pressure indicator	REF 56-02-050	REF 56-02-100	REF 56-02-300	5
With Luer Lock connection	-	-	REF 56-01-300	5

Additional information



Application video

The medical devices in this advertising material are manufactured without the use of natural rubber latex, unless otherwise specified.
The medical devices in this advertising material do not contain phthalates which require labelling according to CLP Regulation (EC) 1272/2008.